

MACHEREY-NAGEL



CHROMABOND®

Promotion 2021

Chromatography



Modern polymeric SPE phases

- CHROMABOND® HLB
for the extraction of medium to polar analytes
from polar matrices
- CHROMABOND® HR-Xperts
for the extraction of hydrophobic and ionic
analytes from polar matrices

Promotion*
Oct. 1st – Dec. 31st,
2021

MACHEREY-NAGEL

www.mn-net.com



CHROMABOND® Promotion 2021

Features of CHROMABOND® HLB

- Applicable for a wide range of analyte polarities
- Enhanced retention for polar compounds
- High loadability and outstanding performance
- Water wettable – even if bed runs dry, SPE can be continued
- Two particle sizes (60 µm and 30 µm) available
- The alternative to Oasis® HLB

Features of CHROMABOND® HR-Xperts

- RP and different mixed-mode polymer phases in two particle sizes (45 µm and 85 µm)
- Broad spectrum of applications with special suitability for enrichment of pharmaceuticals from biological matrices
- High purity
- High loadability and outstanding elution properties

Ordering information

REF	Phase	Description	Volume	Adsorbent weight	Particle size	Pack of	Price per pack
N-vinylpyrrolidone-divinylbenzene (NVP/DVB) resin							
730921	CHROMABOND® HLB	Hydrophilic-lipophilic balance	1 mL	30 mg	60 µm	30	
730921P30	CHROMABOND® HLB	Hydrophilic-lipophilic balance	1 mL	30 mg	30 µm	30	
730923	CHROMABOND® HLB	Hydrophilic-lipophilic balance	3 mL	60 mg	60 µm	30	
730923P30	CHROMABOND® HLB	Hydrophilic-lipophilic balance	3 mL	60 mg	30 µm	30	
730924	CHROMABOND® HLB	Hydrophilic-lipophilic balance	3 mL	200 mg	60 µm	30	
730924P30	CHROMABOND® HLB	Hydrophilic-lipophilic balance	3 mL	200 mg	30 µm	30	
730926	CHROMABOND® HLB	Hydrophilic-lipophilic balance	6 mL	200 mg	60 µm	30	
Polystyrene-divinylbenzene resin							
730936	CHROMABOND® HR-X	Hydrophobic resin	3 mL	60 mg	85 µm	30	
730936P45	CHROMABOND® HR-X	Hydrophobic resin	3 mL	60 mg	45 µm	30	
730931	CHROMABOND® HR-X	Hydrophobic resin	3 mL	200 mg	85 µm	30	
730931P45	CHROMABOND® HR-X	Hydrophobic resin	3 mL	200 mg	45 µm	30	
730939	CHROMABOND® HR-X	Hydrophobic resin	6 mL	500 mg	85 µm	30	
730956P45	CHROMABOND® HR-XC	Strong mixed-mode cation exchanger	3 mL	60 mg	45 µm	30	
730952	CHROMABOND® HR-XC	Strong mixed-mode cation exchanger	3 mL	200 mg	85 µm	30	
730950P45	CHROMABOND® HR-XA	Strong mixed-mode anion exchanger	3 mL	60 mg	45 µm	30	
730951	CHROMABOND® HR-XA	Strong mixed-mode anion exchanger	3 mL	200 mg	85 µm	30	
730735P45	CHROMABOND® HR-XCW	Weak mixed-mode cation exchanger	3 mL	60 mg	45 µm	30	
730739	CHROMABOND® HR-XCW	Weak mixed-mode cation exchanger	3 mL	200 mg	85 µm	30	
730747P45	CHROMABOND® HR-XAW	Weak mixed-mode anion exchanger	3 mL	60 mg	45 µm	30	
730748	CHROMABOND® HR-XAW	Weak mixed-mode anion exchanger	3 mL	200 mg	85 µm	30	

* Special promotion prices only valid for the mentioned products. Offer valid from October 1st – December 31st, 2021

Registered trademarks

Oasis® Waters Corp. (USA)



Labquip (Ireland) Limited
 12 The Business Centre, Fonthill Industrial Park,
 Clondalkin, Dublin 22, D22 X8P5, Ireland
 T: +353 (0)1 643 4586
 E: labquip@labquip.ie
 Website: www.labquip.ie

www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG
 Valencienner Str. 11
 52355 Düren · Germany

DE Tel.: +49 24 21 969-0 info@mn-net.com
 CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com
 FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com
 US Tel.: +1 888 321 62 24 sales-us@mn-net.com